**VOCO**

# Structur Premium/QM

**(GB)**

Instructions for use

**CE** 0482

**Structur Premium/QM** is a cold polymerizing paste-paste system for the production of temporary crowns, bridges, inlays and onlays. **Structur Premium/QM** has excellent mechanical properties and is fluorescent. **Structur Premium** is available in 7 shades: A1, A2, A3, A3.5, B1, B3 and BL. **Structur Premium QM** is available in the shade I (incisal).

**Directions:****Structur Premium in cartridges:**

Insert **Structur Premium** cartridge into the **VOCO** Dispenser (type 2) or a compatible extruder gun with a suitable piston. Due to technical reasons, slight differences in the filling level between the cartridge chambers may occur. In this case, remove cartridge cap and squeeze the material out until the material is evenly delivered from both chambers. Then insert a mixing tip type 6 and lock with a 90° turn clockwise. By pressing the lever, the material is extruded and automatically mixed in the correct ratio.

Do not remove the mixing tip from the cartridge after use. Replace the used mixing tip with a new one just before the next application. When changing the mixing tip, check the opening of the cartridge to confirm that the material can flow freely.

**Structur Premium QM in the QuickMix syringe:**

Remove the cap from the syringe after turning it counter-clockwise. Attach a mixing tip type 10. Match the nocks on the syringe to the mixing tip and secure the mixing tip with a 90° turn clockwise. The material is automatically mixed in the correct ratio when dispensing and can be applied directly.

All QM syringes come with a detachable piston.

The piston can simply be pulled from the QM syringe and replaced without causing any material leakage. The plugs in the syringe chambers will stay in place, preventing any cross-contamination of the material.

**Advice on application:**

Level out any undercut in the preparation. Do not use cements containing eugenol for this purpose.

Prepare an impression if prefabricated crowns or similar restoratives are not being used.

Prepare the impression by levelling undercuts and cutting drainage channels.

Enlarge the impression if the wall thickness of the temporary restoration is too thin, taking care to remove the septums in the interdental embrasures. Clean and air dry the impression afterwards.


Apply **Structur Premium/QM** directly into the impression. Do not overfill.

The impression must be applied in the mouth within 30 s. Remove the temporary from the mouth during the elastic phase, 1:00 min - 1:30 min after mixing.

The setting process can also be monitored by testing the excess material in the mouth.

The setting time of non-used material (in mixing tips) **cannot** be compared to the setting time of the material in the mouth.

At a room temperature of 23°C and 50% humidity, working times are as follows:

Slight differences may result from different room temperatures.

After ca. 4 min, take the temporary out of the impression and wipe dry. If necessary, a solvent such as alcohol may be used. This removes the thin inhibition layer usually caused by oxygen and allows easier finishing. The surface of the temporary restoration should appear silky-matt.

Carefully trim any excess material and roughen areas near the gingiva. Cement the temporary restoration using an eugenol-free temporary cement.

**Making temporaries with the 2 layer technique:**

After the mould has been taken and prepared, apply **Structur Premium QM I** directly onto the bottom of the mould and smooth it out toward the cervical area. Apply the selected shade of **Structur Premium**. No additional waiting time is required.

**Warnings and precautions:**

- **Structur Premium/QM** contains methacrylates, amines, terpenes, benzoyl-peroxide, BHT and should not be applied in cases of known allergies to these ingredients.

- Before making the temporary, prepared cores should be insulated to promote desensitization (fluoride varnishes, core varnishes, etc.).

- In order to avoid fractures or injuries from debris, the temporary should not be subjected to excessive loads (molar area).

- Linings and cores built up with resin-based materials (composites) should be sufficiently insulated to prevent removal with the temporary (e.g. glycerine).

- Inadequate hygiene or frequent use of certain mouthrinses may lead to discolouration.

- Finishing should only be started after material has completely set (approx. 4 min after mixing).

- Restoration of defective temporaries:

If the temporary breaks despite the extreme stability of **Structur Premium/QM**, the following procedures are recommended:

- **Air bubbles:** fill with **Structur Premium/QM**.
- **In case of fracture immediately after manufacturing the temporary:** re-connect the fracture with freshly mixed **Structur Premium/QM**.
- **Repairing fractured temporaries that have been worn:**

Roughen the fracture with a bur and place undercuts. Apply bond on the prepared fractures and cure according to the respective instructions for use. Apply a flowable composite on the prepared areas, assemble the fractured pieces and cure according to the instructions for use.

**Storage:**

Store at temperatures between 4°C - 23°C.

**Protect from bright light and store in closed packages, drawers, etc.**

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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